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Instructions for use

Patient target group

All persons who need individual trays are required for treatment.

Intended user

Dentista

Composition Qualitative data:

Oligomero: urethane acrylate, vinyl ester resin, 2,2-dimethoxy-1,2-diphenylethan-1-one or by 2,2-dimethoxy-2-phenylacetopheno, (+)-2,3-dioxo-1,7,7-trimethyl-bis(cyclo-(2,2,1)-heptano di-borano-2,3 dione, difenil(2,4,6-trimetilbenzolo)fisofina, 2-hidrox-2-metil-1-phenyl-prop-1-one, 2-ettesil-4-(dimetilammino)benzoato, 2,2 Dietossacetofeno, vetro sodocalcico, silice amorfa, coloranti (Pigment Red 146 = CI 12485 / Pigment Blue 29 = CI 77007 / D&amp;C Red 30 = CI 73360 / Pigment Blue CI 74160 / Pigment Violet 19 = CI 73900), flavorings (Optamint).

Quantitative data:

In 100 g there are 19 g acrylates/methacrylates and 80 g of vetro.

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Elite LC Tray / Elite LC Tray Round polymerises with light to form a solid and compact polymer.

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Avoid any contact between the eyes and skin and the uncured material, i.e. also with the inhibited layer.

Always use a diatherm when trimming the tray!

Preparation/Manufacture

- Mark the outline of the future tray on the model.

- Block out the undercuts on the model.

- To create a spacer for the impression material, block out the entire model (e.g. with a sheet of wax) and cut out stops.

- If the plaster is very dry, soak the model briefly in water or apply dye isolation or an alginate isolation.

- Remove the baseplate from the box and close it.

- Once the model has been blocked out, adapt the baseplate to it and cut to shape using, e.g. a scalpel.

- Use the remaining material to form a handle and attach it to the tray.

- Place the model and tray in a light-curing unit and cure without vacuum.

- Remove the tray from the model and polymerize once again from the underside.

- If the tray has been cured without a vacuum, use alcohol to remove the smear layer.

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The top of the tray, on the model 3-5 minutes

The underside of the tray, off the model 3-5 minutes

Trimming:

- Trim the custom tray with tungsten carbide burs and smooth it with sandpaper.

Hygiene:

Immerse impressions taken with Elite LC Tray / Elite LC Tray Round customs trays in a disinfectant solution based on glutaraldehyde. The time depends on the disinfectant manufacturer's directions.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

Caution

Adverse effects are not known.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

The product should be stored between 0°C - 22°C

Technical data

Flexural strength min. 130 N/mm<sup>2</sup>

Minimum shelf life

Do not use Elite LC Tray / Elite LC Tray Round after the expiry date.

Variants

- Elite LC Tray [REF] D500100 White [REF] D500101 Pink [REF] D500102 Blue

- Elite LC Tray Round [REF] D500052 Pink

Disposal

Uncured material can be polymerised with light (sunlight or artificial light) and then disposed of as plastic waste.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

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